

Agilent Scan Control V9.1 Software Installation Notes

To install Scan Control software, double click on Setup.exe:

- You must be logged in as an Administrator to perform this installation.
- Setup.exe will first install Agilent Installation Qualification Tool (IQT).
- It will then install Agilent Microarray Scan Control software.
- Agilent recommends performing installation to the default folder.

Notes:

- You do not need to uninstall previous versions of the software.
- If you need to reinstall the current version, select the Repair option when prompted. This will perform a full reinstall.
- Configuration data, such as Scan Regions and Scan Protocols, will be preserved during upgrades and reinstallations. Instrument calibration parameters are stored on the instrument itself, and are not affected by reinstallation.
- You will be prompted several times by Windows User Access Control to permit installation to proceed. Please click Yes when prompted.
- When you first connect to the scanner after installing a software update, you will be prompted to install firmware on the instrument. Please follow the instructions presented; you will need to restart the instrument and Scan Control after firmware is updated.

Restoration of Firewall Settings

Correct operation of the scanner requires that Scan Control be allowed full access through the Windows Firewall. This rule is set by the installer. If the firewall rules are reset for any reason, you can re-apply the Scan Control rules by double clicking on SetFirewallRules.bat in the scanner executable directory (usually C:\Program Files (x86)\Agilent\ScanControl).